


CHEMILUMINESCENT KIT for HRP DETECTION - WESTERN BLOTTING

PRODUKT	KATALOG. ČÍSLO	VELIKOST A BALENÍ		CENA [CZK]	DETEKCE JIŽ OD
WESTAR SUN	XLS063,0250	250 mL kit	1 Kit	2 671 Kč	mid-picograms
WESTAR NOVA 2,0	XLS071,0250	250 mL kit	1 Kit	3 075 Kč	low-picograms
WESTAR ANTARES	XLS142,0250	250 mL kit	1 Kit	4 018 Kč	mid-femtograms
WESTAR ETAC ULTRA 2.0	XLS075,0100	100 mL kit	1 Kit	6 983 Kč	mid/low-femtograms
WESTAR SUPERNOVA	XLS3,0100	100 mL kit	1 Kit	8 330 Kč	low-femtograms
WESTAR HYPERNOVA	XLS149,0020	20 mL	1 Kit	3 883 Kč	attograms



WESTAR SUPERNOVA - Cyanagen

Clarity Max™ - Bio-Rad

Amersham™ ECL Select™ - GE Healthcare

SuperSignal™ West Femto - Thermo Scientific™

Western blotting detection of HDAC-1 on Hela cell lysate with WESTAR SUPERNOVA and its competitors

Sample: Two-fold dilution series of Hela whole cell lysate (abcam®) from 2.5 µg to 0,039 µg of total protein

Membrane: Trans-Blot® Turbo™ Mini Nitrocellulose Transfer Packs (Biorad)


Blocking: 2% ECL™ Blocking Agent (GE Healthcare) in PBS-T

Primary antibody: Rabbit-anti Human HDAC-1 (abcam®) 1:15000

Secondary antibody: Goat anti-Rabbit IgG HRP (2mg/mL) (abcam®) 1:300000

Imaging: ImageQuant™ LAS 4000 (GE Healthcare)

Exposure time: 120 seconds.



WESTAR SUPERNOVA Signal Duration Comparison with Competitors

Comparison of signal intensities at time points up to 20 hours post substrate addition.

Exposure time is 120 seconds for each time point (0-1-2-5-8-11 hours).

CHEMILUMINESCENT KIT for HRP DETECTION - WESTERN BLOTTING - SINGLE SOLUTION

PRODUKT	KATALOG. ČÍSLO	BALENÍ	CENA [CZK]	DETEKCE JIŽ OD
WESTAR-ONE BASIC	XLSU177,0250	250 mL kit	2 805 Kč	mid-picograms
WESTAR-ONE PLUS	XLSU178,0250	250 mL kit	3 479 Kč	mid-femtogram
WESTAR-ONE EXTREME	XLSU180,0125	125 mL kit	8 330 Kč	low-femtogram

ECL in TABLETS for WESTERN BLOTTING

PRODUKT	KATALOG. ČÍSLO	BALENÍ	CENA [CZK]	DETEKCE JIŽ OD
WESTAR-TAB SUPER	XLSC310,020	20 (equiv to 100ml)	8 330 Kč	low-femtogram

WESTERN BLOTTING - STRIPPING BUFFER SOLUTION

PRODUKT	KATALOG. ČÍSLO	BALENÍ	CENA [CZK]	
RENEW	SBS069,0050	50 mL	919 Kč	
	SBS069,0250	250 mL	2 132 Kč	
	SBS069,0500	500 mL	3 344 Kč	

BCA PROTEIN ASSAY						
PRODUKT	KATALOG. ČÍSLO	VELIKOST A BALENÍ		CENA [CZK]	PRACOVNÍ ROZSAH	OBSAH DOSTAČUJÍCÍ PRO
QPRO - BCA Kit Standard 20 - 2000 µg/mL	PRTD1,0500	Reagent A 500 mL Reagent B 15 mL	1 kit	3 344 Kč	20÷2000 µg/ml	250 test tube or 2,500 microplate assays
µQPRO - BCA Kit Micro 0,5 ÷ 20 µg/mL	PRTD2,0500	Reagent A 240 mL Reagent B 240 mL Reagent C 12 mL	1 kit	6 713 Kč	0,5÷20 µg/ml	480 test tube or 3,200 microplate assays

STAR BEADS UNIVERSAL DNA/RNA EXTRACTION Kit					
PRODUKT	KATALOG. ČÍSLO	FORMÁT	BALENÍ	CENA [CZK]	
STAR BEADS UNIVERSAL DNA/RNA EXTRACTION Kit	SBK270,1X96	bottles	1x96	6 039 Kč	
STAR BEADS UNIVERSAL DNA/RNA EXTRACTION Kit PreFilled Plate	SBK271,1X96PFI	pre-filled plates	1x96	9 543 Kč	
STAR BEADS UNIVERSAL DNA/RNA EXTRACTION Kit PreFilled Plate	SBK273,2x32PFI	pre-filled plates	2x32	6 578 Kč	

STAR BEADS DNA EXTRACTION Kit					
PRODUKT	KATALOG. ČÍSLO	FORMÁT	BALENÍ	CENA [CZK]	
STAR BEADS DNA EXTRACTION Kit	SBK307,1X96	bottles	1x96	6 039 Kč	
STAR BEADS DNA EXTRACTION Kit PreFilled Plate	SBK307,1X96PFI	pre-filled plates	1x96	10 082 Kč	
STAR BEADS DNA EXTRACTION Kit PreFilled Plate	SBK307,2X32PFI	pre-filled plates	2x32	6 848 Kč	

STAR BEADS RNA EXTRACTION Kit					
PRODUKT	KATALOG. ČÍSLO	FORMÁT	BALENÍ	CENA [CZK]	
STAR BEADS RNA EXTRACTION Kit	SBK308,1X96	bottles	1x96	6 039 Kč	



STAR BEADS cfDNA/RNA EXTRACTION Kit					
PRODUKT	KATALOG. ČÍSLO	FORMÁT	BALENÍ	CENA [CZK]	
STAR BEADS cfDNA/RNA EXTRACTION Kit	SBK309,1X96	bottles	1x96	11 564 Kč	

STAR BEADS LYSIS BUFFER					
PRODUKT	KATALOG. ČÍSLO	FORMÁT	BALENÍ	CENA [CZK]	
STAR BEADS Lysis Buffer	SBPLB254,0030	bottle	30 mL	2 132 Kč	
STAR BEADS Lysis Buffer	SBPLB254,0125	bottle	125 mL	6 444 Kč	
STAR BEADS Tissue Lysis Buffer	SBLB309,0020	bottle	20 mL	2 132 Kč	

STAR BEADS SILANOL					
PRODUKT	KATALOG. ČÍSLO	JÁDRO	BALENÍ	CENA [CZK]	
STAR BEADS Silanol	SBB306,6ML	silica	6 mL	4 288 Kč	
STAR BEADS Silanol	SBB306,50ML	silica	50 mL	20 458 Kč	

STAR BEADS PROTEIN G					
PRODUKT	KATALOG. ČÍSLO	JÁDRO	BALENÍ	CENA [CZK]	
STAR BEADS PROTEIN G	SBBG321,5ML	silica	5 mL	15 068 Kč	
STAR BEADS PROTEIN G	SBBG321,50ML	silica	50 mL	99 960 Kč	

Rádi Vám sdělíme na požádání ceny dalších produktů firmy Cyanagen pro

WESTERN BLOTTING, ELISA, CLIA, FIA, FLUORESCENT LABELING, FLUORESCENT GEL STAINING, IMAGING, DNA/RNA/PROTEIN ANALYSIS, FLOW CYTOMETRY atd.

Biologicals s.r.o.


U Trati 44a, Praha 10, 100 00

E-mail: info@biologicals.cz

Tel: +420-274776731

E-shop: www.biologicals.cz

Webové stránky výrobce: <https://www.cyanagen.com/>



IR-BLOT Secondary Antibodies